














Find your pipette model and size in chart and appropriate tip size at the top of the column

Pipette Model	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL	5 mL Tips	10 mL Tips
Ovation ESC, E8 & E12	0.5 - 20 µL 	2 - 125 µL  5 - 250 µL 	25 - 850 µL  25 - 1250 µL 		
Ovation QS & M	0.2 - 2 µL  1 - 10 µL  2 - 20 µL 	10 - 100 µL  20 - 200 µL 	100 - 1000 µL 		
Ovation F1 & F2	1 to 20 µL	21 to 200 µL	201 to 1000 µL		
Ovation Macro5 & Macro10				100 µL - 5 mL 	200 µL - 10 mL 
All MLA models	n/a	2 to 200µL	201 to 1000 µL	1 mL to 5 mL	6 mL to 10 mL



Description	Product #		
	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL
<b>VistaRak™</b> 192 tips/rack, 5 racks	4060-1002	4060-2004	4060-3004
<b>VistaRak, Low Retention</b> 192 tips/rack, 5 racks	4070-1002LR	4070-2004LR	4070-3004LR
<b>VistaStak™</b> 192 tips/layer 3 layers for Micro (25 µL) 5 layers for Small (250 µL) tips 3 layers for Large (1250 µL) tips	4060-9024	4060-9025	4060-9026
<b>Stacked Rack</b> 200 tips/layer 5 layers for Small (250 µL) tips 3 layers for Large (1250 µL) tips		9025 <sup>1</sup>	9026 <sup>1</sup>
<b>VistaBulk™</b> 1000 tips/bag	4058-1000	4058-2000	4058-3000
<b>Protectainer™ Bulk Pack</b> 1000 Small (250 µL) tips 750 Large (1250 µL) tips		4025	4026

1. Not for use with Ovation Multichannel pipettes.

Non-Sterile Tips

## Sterile

Description	Product #		
	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL
<b>VistaRak, Sterile</b> Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4060-1032	4060-2132	4060-3132
<b>VistaRak, Sterile , Low Retention</b> Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4070-1032LR	4070-2132LR	4070-3132LR
<b>VistaTip, Individually wrapped, Sterile</b> Certified RNase/DNase free, non-pyrogenic 200 tips		2027	2028



VistaTip

## Sterile Filtered

Description	Product #		
	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL
<b>VistaRak, Sterile, Filtered</b> Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4060-1332 <sup>1</sup> 4060-1333 <sup>2</sup>	4060-2332	4060-3332
<b>VistaRak, Sterile , Filtered, Low Retention</b> Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4070-1332LR <sup>1</sup> 4070-1333LR <sup>2</sup>	4070-2332LR	4070-3332LR

1. Not for use with Ovation QS and Ovation M 0.2-2 µL  
2. 2 µL tips for use with Ovation QS and Ovation M 0.2-2 µL



VistaTip (Macro tip)

## Macro Volume Tips

Description	Product #	
	Macro Tips 5 mL*	Macro Tips 10 mL
<b>Box, Non-Sterile</b> 100 tips/box for 5 mL 35 tips/box for 10 mL	9048	9050
<b>VistaClear™ Box, Non-Sterile</b> 35 tips/box		4058-6100
<b>VistaBulk, Non-Sterile</b> 250 tips/bag for 5 mL 100 tips/bag for 10 mL	4058-5000	4058-6000
<b>VistaClear Box, Sterile</b> Certified RNase/DNase free, non-pyrogenic 60 tips/box for 5 mL 35 tips/box for 10 mL	4058-5102	4058-6102
<b>VistaClear Box, Sterile, Filtered</b> Certified RNase/DNase free, non-pyrogenic 60 tips/box for 5 mL 35 tips/box for 10 mL	4058-5133	4058-6133
<b>VistaTip, Individually wrapped, Sterile, Filtered</b> Certified RNase/DNase free, non-pyrogenic 50 tips	4058-5332	4058-6332

\*5 mL macro tips are graduated at 0.1 mL increments.

VistaClear Box  
(Macro tips)Box  
(Macro tips)